



University of Huddersfield Repository

Bills, Paul J.

Introduction to Computed Tomography as an aid to Forensic Analysis

Original Citation

Bills, Paul J. (2016) Introduction to Computed Tomography as an aid to Forensic Analysis. In: Age Estimation Workshop, 13th - 14th May 2016, University of Huddersfield, UK. (Submitted)

This version is available at <http://eprints.hud.ac.uk/id/eprint/28556/>

The University Repository is a digital collection of the research output of the University, available on Open Access. Copyright and Moral Rights for the items on this site are retained by the individual author and/or other copyright owners. Users may access full items free of charge; copies of full text items generally can be reproduced, displayed or performed and given to third parties in any format or medium for personal research or study, educational or not-for-profit purposes without prior permission or charge, provided:

- The authors, title and full bibliographic details is credited in any copy;
- A hyperlink and/or URL is included for the original metadata page; and
- The content is not changed in any way.

For more information, including our policy and submission procedure, please contact the Repository Team at: E.mailbox@hud.ac.uk.

<http://eprints.hud.ac.uk/>

Introduction to Computed Tomography as an aid to Forensic Analysis

Dr Paul Bills

EPSRC Centre for Innovative Manufacturing in Advanced Metrology

Centre for Precision Technologies

University of Huddersfield, United Kingdom

Agestimation Workshop, Huddersfield 13-14 May 2016

EPSRC Centre for Innovative Manufacturing in Advanced Metrology.

Nationally funded, centre of excellence in advance metrology. Based at the University of Huddersfield's Centre for Precision Technologies, with an international reputation in precision engineering, metrology research and standards development.


Key areas of research are:

- Additive Manufacturing
- Software Development
- Hardware Applications
- Surface Measurement & Applications
- Ultra Precision Manufacturing
- Industrial Metrology



What are X-rays?


- X-rays are electromagnetic radiation just like visible light, infra-red light, ultra-violet light and radio waves, but with a much shorter wavelength than any of these.




Typical μ CT X-ray sources produce energies in the range 30-450keV (in red).

How do we generate X-rays?

- We generate X-rays by firing electrons at high speed on to a metal target.
- Electrons are produced from a hot filament (like a light bulb).
- They are accelerated using a high voltage into a beam tube.
- They travel at up to 80% the speed of light (giving them energies of 30 - 450keV).
- They are focused by a magnetic lens into a small spot (1 – 5 μ m) onto a metal target.
- The sudden deceleration of the charged electrons when they hit the target produces 99.3% heat and 0.7% X-rays.




X-ray Source: Reflection target




Note: Reflection target limits the maximum magnification as the target is behind the vacuum window

How do we get an X-ray image?



X-rays travel in straight lines and pass right through the sample. However, some of the X-rays are absorbed by the sample and so the intensity of the X-rays is reduced forming a shadow image.

How do we get a magnified image?




Just like light, X-rays travel in straight lines.

Unlike light, we cannot use a lens, so we use geometric magnification.


The magnification is increased by moving the sample closer to the X-ray source (and vice versa).

What is Computed Tomography?

- Computed Tomography (or CT) is the process of imaging an object from all directions using penetrating radiation (e.g. X-rays) and using a computer to reconstruct the internal 3-D structure of the object from the intensity values in the projected images.
- It is the process used in a medical CT scanner, though in our case we keep the source and detector stationary and rotate the object. Hospital patients might complain if we did this to them!




How does CT actually work?




How does CT actually work?

CT requires us to penetrate the object with X-rays from all directions:




How does CT actually work?

CT requires us to penetrate the object with X-rays from all directions:








How does CT actually work?

CT requires us to penetrate the object with X-rays from all directions:




How does CT actually work?



From thousands of images like these a computer algorithm generates a 3D volume which can be sliced in software to reveal the internal structure of the object.

Changing pixels into voxels


Projecting the volume elements (voxels) onto the picture elements (pixels)




Each voxel is projected on to the pixel that a line from the source through that voxel hits. The intensity of this pixel is added into the voxel. This is repeated for all images. So for 3000 images, each of the billion or so voxels is processed 3000 times. That's why it can take a long time to reconstruct high-resolution CT volumes!

Hip replacement adhesion







Bone degradation modelling



Cremation offering analysis



Is rice nice at twice the price?



Subtitle: The evils of weevils

Fly anatomy






Fly anatomy



Fly anatomy



CT for ballistic analysis



The authors gratefully acknowledge the European Research Council (ERC-ADG-228117) funding and UK's Engineering and Physical Sciences Research Council (EPSRC) funding of the EPSRC Centre for Innovative Manufacturing in Advanced Metrology (Grant Ref: EP/I033424/1).



EPSRC


Engineering and Physical Sciences
Research Council

Acknowledgements

- Katie Addinall
- Chris Dawson



Any Questions?



Contact:
Paul Bills
p.j.bills@hud.ac.uk
+441484472769